



# SureDental™

*Presents  
Proven Products from*

**parkell**  
AUTHORIZED DEALER *inc.*

## Blu-Mousse *Bite Registration*

Exceptional accuracy.  
30-second set  
Plaster-like hardness (90-durometer).  
Also available in peppermint scent. (S448S)  
Dependable, no-slump consistency.  
Does it all from conventional bite registration techniques to modern-day digital impressions.

Specify Setting Time

Blu-Mousse Classic Regular Set (2min Set)

Blu-Mousse Classic Super Fast Set (30second Set)

Blu-Mousse 60 (60second Set)

## Brush & Bond

*Dentin/Enamel Bonding Liquid*



**\$127.95**  
3 ml Btle

Bonds direct and indirect composite restorations to all tooth surfaces.  
Permanently seals the tubules to prevent sensitivity.  
Protects the pulp against microleakage and bacterial invasion.  
Creates an acid and wear-resistant shield for the tooth.  
Doubles the surface hardness of the preparation.  
Cures with all currently available curing lights.  
Low film thickness won't affect crown fit.  
Also available— Brush&Bond® UNIVERSAL  
Brush & Bond Kit

One 3ml bottle of Dentin/Enamel  
Bonding Liquid  
100 Standard Activator Brushes  
Mixing well

**\$217.95**  
Complete Kit

## Cinch

**Polyvinyl Impression Material**

Exceptionally accurate, hydrophilic and strong.  
Use as a final impression material or an alginate substitute.  
4 min. set, 3 min. working time.  
Cinch Platinum is also available in tubes and a 380ml super cartridge  
Cinch is available in a wide variety of viscosities, final hardness and setting times.



# SureDental

The Gold Standard  
in Bite Registration



**Buy 3  
Get 1  
FREE**


**\$64.95**  
Each

Two 50ml split cartridges (100ml ttl.) Mixing tips.

**Net  
\$48.71**



Blu-Mousse Blu-Mousse 60 Green-Mousse



**Buy 3  
Get 1  
FREE**

**\$84.95**  
4 50ml

1.888.277.2335



**The TurboVue®** features a light source built into the Complete Antimicrobial Delivery System handpiece allowing a significant amount of light to emit through the 30K, light-transmitting ultrasonic inserts

### Integra (Self contained Ultrasonic Scaling System)

Features a stable reservoir bottle design allowing the operator to gently move the scaler from operator to operator (on a mobile cart) without spilling the liquids from the bottles.



### TurboSENSOR (Ultrasonic Scaler)

**THE WORKHORSE** – The Hardest Working and Most Practical Ultrasonic Scaler. Parkell's popular TurboSENSOR® Ultrasonic Scaler combines true versatility with intuitive operation. Its smart-circuitry senses between a 25KHz or 30KHz tip and automatically switches operating frequency to match the insert in the handpiece.

**Buy a TurboSensor®, Get a Free Insert\***

## Materials



**Stops gingival bleeding and seepage that may interfere with impression taking**  
Great for tissue management when seating restorations, placing rubber dam clamps (controlling bleeding), bleaching teeth (to manage moisture) and restoring subgingival cavities  
Fast-acting, Fraction of the cost of all other retraction materials  
Light green color makes the material easy to detect  
Also available in Unit Dose Capsules

Photo Credit: Henry Long, DMD

*We've used Dryz for a few years and I'm always impressed at how it stops bleeding almost instantly. I usually use laser retraction but when there is bleeding, I will often Dryz first. When washed off, I rarely have to pick up the laser. The stuff is just plain amazing. —Steven H. Schwartz, D.M.D.*

### C & B Metabond Quick Luting Cement

Cement crowns on short or overly tapered abutments. Re-cement bridges or posts that repeatedly loosen due to compromised retention or patient habits. Perform fast, chairside repairs such as mending broken abutments, repairing porcelain fractures, splinting periodontally compromised teeth, adding teeth to removable prostheses, or bonding directly to composite or acrylic-based prostheses. Place direct pulp-caps, repair root perforations in "hopeless cases" and so much more...



**CALL**

- DCG30Uni-PF
- DCG30Uni
- DLH30Tri-PF
- DLH30Tri
- DUI30
- DUI25
- DTI30
- DTI25
- DPR30
- DPR25
- DPL30
- DPL25
- DPI30
- DPI25
- DBI30
- DBI25
- DCG30Tri-PF
- DCG30Tri
- DCG30Slim-PF
- DCG30Slim
- DSG30U-P
- DSG30U
- DSG30T
- DSG30R-P
- DSG30R
- DSG30P
- DSG30P-P
- DSG30L-P
- DSG30L
- DLH30T
- DLH30U
- DLH30P

**CALL**

## Materials



**NEW**

### AccuFilm® Fastcheck— Buy 3 Get 1 Free.\*

Built-in holders eliminate the need for forceps. Consistently produces sharp, easy-to-interpret marks. Simultaneous marking of both arches with its double-sided film. Easily identify high contacts making for easier adjustments. Incredibly thin film for accurate exposure of high occlusal contacts. Works for a single-tooth or entire quadrant. For use on wet or dry surfaces.

**Net \$26.21**

**\$34.95 ea**

### Blu-Mousse®— Buy 3 Get 1 Free.\*

Blu-Mousse® & Green Mousse® Bite registrations, matrices, model bases, bite plates and now CAD scanning capabilities—the uses for Blu-Mousse® just keep growing! Exceptional accuracy. 30 Second, 70-second & 2 Minutes set.

**Net \$48.71**

**\$64.95 ea**

### Cinch®— Buy 3 Get 1 Free.\*

The unique properties are exceptionally adaptable and capable of capturing all of the finer details of a preparation area

**Available in Light, Medium, Heavy and Auto mix putty**

**Net \$63.71**

**\$84.95 ea**

### Dryz®— Buy 3 Get 1 Free.\*

Dryz stops gingival bleeding and seepage that may interfere with impression taking. Great for tissue management when seating restorations, placing rubber dam clamps (controlling bleeding), bleaching teeth (to manage moisture) and restoring subgingival cavities

**HemoStatic Retraction Paste**

**7 x .5ml Syringes + Tips \$104.95 ea**

### Brush&Bond®—

Tried & true—the easy-to-use, high-strength bonding agent that protects and desensitizes! Use it to bond direct and indirect restorations to all tooth surfaces, seal crown preps, assist with core build-ups, and act as a relief agent for patients with dentin hypersensitivity

**Liquid Only 3ml 127.95 Complete Kit \$217.95 ea**



**SureDental™**

**1.888.277.2335**