

## Endo Sealer Solvent<sup>TM</sup>

The first endodontic irrigant designed for bioceramic sealer retreatment and calcium hydroxide removal.

- **Targeted Efficacy:** Engineered to selectively remove sealer and calcium hydroxide while being gentle on dentin
- **Versatile:** Suitable for all bioceramic sealers and calcium hydroxide
- **Time-Efficient:** Streamlines two visit endodontic procedures & retreatment procedures

Endo Sealer Solvent<sup>TM</sup> Bottle

317033 | 1 x 8oz (240mL) Bottle

Buy 1 Kit  
Get:

(1) Appli-Vac<sup>TM</sup> 1" 30GA  
Irrigating Tips (100pk)  
315530



## EndoUltra<sup>®</sup> Ultrasonic Activator

Cordless Ultrasonic Irrigation

- Better penetration of irrigants throughout the canal
- Better acoustic streaming with cavitation
- Better cleansing in less time

Buy 1 Kit, Get  
(1) EndoUltra<sup>®</sup>  
Activator Tips  
(3pk) 407706  
**FREE**



EndoUltra<sup>®</sup> Cordless Ultrasonic Activator

- 407702
- 1 x Cordless handpiece
  - 3 x Titanium 20/02 Activator Tips
  - 3 x Autoclavable Handpiece Sleeves
  - 100 x Protective Barriers
  - 1 x Tip Wrench
  - 1 x Charger

## Appli-Vac<sup>TM</sup> 1" Irrigating Tip

- Luer Lock Hub
- Slotted and side-vented for safe irrigation

Appli-Vac<sup>TM</sup> 1" Irrigating Tip

- 315523 Appli-Vac<sup>TM</sup> 1" Irrigating Tip 23ga 100pk  
100 x 23ga (Blue)
- 315527 Appli-Vac<sup>TM</sup> 1" Irrigating Tip 27ga 100pk  
100 x 27ga (Yellow)
- 315530 Appli-Vac<sup>TM</sup> 1" Irrigating Tip 23ga 100pk  
100 x 30ga (Purple)

Buy 3  
Get 1  
**FREE**



# Q2 Canada Offers 2026

Free goods will ship with order | Offers run 4-1-2026 through 6-30-2026

Use code: **VAQ2\_ENDO\_C**



**VistaApex**  
SOLUTIONS  
Value. Added.

1.888.277.2335

## SmearOFF™

### 2-in-1 EDTA/CHX Mix

- EDTA / CHX mix formulated with wetting agents and surface modifiers
- Will NOT form a precipitate when used with NaOCl
- Does not require a rinse



### SmearOFF™ 16oz Bottle

317111 | • 1 x 16oz Bottle

Buy 1  
Get

Chlor-XTRA™  
16oz Bottle  
503800



## EDTA

### 17% Solution

Specifically engineered with LOWER pH for optimal performance! For best outcomes, a LOWER pH is more important than concentration

- Opens dentin tubules
- Removes smear layer and dentin mud
- Prepares dentin walls for better ADHESION of filling materials



### EDTA 17% Solution 16oz Bottle

317011 | • 1 x 16oz Bottle

Buy 1  
Get

6% Sodium  
Hypochlorite  
16oz Bottle  
317008



## CHX-Plus™ Solution

- 2% Chlorhexidine Gluconate (CHX) enhanced with powerful wetting agents and proprietary surface modifiers
- Added surfactants lower surface tension, improving depth of penetration
- Promotes adhesion of sealers and obturation materials



### CHX-Plus™ Solution 16oz Bottle

503900 | • 1 x 16oz Bottle

Buy 1  
Get 1  
FREE

**VistaApex**  
SOLUTIONS  
Value. Added.

2200 South St, Racine, WI 53404  
262-636-9755

VistaApex Solutions Promotional Terms & Conditions

Offer(s) valid only in Canada. Participation is not permitted by customers enrolled in any other pricing program, including but not limited to government entities, group purchasing organizations (GPOs), Dental Service Organizations (DSOs), group practices, special markets, community health centers, schools, institutions, or contracted discount programs. Offer may not be combined with any other promotion, rebate, discount, or special pricing agreement. Proof of purchase may be used for internal reporting and marketing purposes. All company and product names are trademarks of VistaApex Solutions, its affiliates, or related companies unless otherwise noted. Trademarks may not be registered in all jurisdictions. Purchases must be made through an authorized VistaApex Solutions distributor between April 1–June 30, 2026.